

Boil Out Unit Wapo-Ex 8 III



USER MANUAL

Dear customer,

Thank you for choosing a product from the Wassermann range. Wassermann Dental-Maschinen incorporates the highest standards of quality and the latest technology. In order to enjoy maximum performance and years of trouble-free operation, please read this user manual carefully before you connect this device and start work, and operate the device according to the recommended guidelines. The operation safety and the functionality of this device can only be guaranteed if you follow both the general safety guidelines and the applying laws to prevent accidents as well as the precautions given in this user manual. We are not liable for any damages which occur due to inappropriate usage or faulty operation of this device.



Make sure that anyone using this device has read and understood this user manual.

Keep this user manual in a safe place where it can be referred to as required at any time.



The unit complies with the relevant EU guidelines.



The unit is subject to the EU guidelines 2012/19/EU (WEEE Directive).

Company address:

Wassermann Dental-Maschinen GmbH

Rudorffweg 15-17

21031 Hamburg, Germany

Phone : +49 (0)40 730 926 -20/ -24 Fax.: +49 (0)40 730 37 24

export@wassermann-dental.com

www.wassermann.hamburg

Contents

1 Features.....	4
2 Safety Symbols used in this Manual.....	5
3 Safety Guidelines.....	5
4 Responsibility for Operation or Damage.....	6
5 Application.....	6
6 Before Starting.....	6
6.1 Transport.....	6
6.2 Installation.....	7
6.3 Storage.....	7
7 Installation/ Commissioning/ Menu Functions.....	8
7.1 The Set-Up Menu.....	9
7.2 The Set-Up Menu Page 1.....	10
7.3 The Set-Up Menu Page 2.....	10
7.4 Description of Functions 5 - 14.....	11
8 Operation.....	15
8.1 General Operating Instructions.....	15
8.2 Operation.....	16
8.3 Hand Shower.....	16
9 Troubleshooting.....	18
10 Maintenance.....	19
10.1 Cleaning.....	19
10.2 Maintenance.....	20
10.3 Repairs.....	20
10.4 Spare Parts.....	20
10.5 Service Hotline 0049 (0)40 730 926 -20/ -24.....	20
10.6 Scope of Delivery/ Accessories.....	20
10.7 Warranty.....	21
11 Technical Data.....	21
12 Disposing of the Unit.....	22
12.1 Information on Disposal for Countries within the EU.....	22
13 EU Declaration of Conformity.....	23

1 Features

The approved wax-scalding unit Wapo-Ex 8 III of the WASSERMANN-brand - nowadays equipped with an easy to clean multi-function display - is a very compact table unit designed for wax-scalding of up to 8 flask halves. Its immersed full-conical spray system simplifies handling and prevents accidental erosion. Its longevity is supported by the silent, joint-free bronze pump and a housing and tank made of stainless steel.

To protect the heating element the wax-scalding unit is equipped with a water-shortage protection mechanism. Furthermore it has a low-lying perforated screen cover as well as a hand shower made of non-heat-conducting plastic for manual scalding of small parts.

- compact tabletop unit for scalding up to 8 flask halves
- informative multifunctional touch display for easy operation
- integrated timer
- thermostat-protected, gasketless immersion pump in bronze; maintenance-free
- energy-saving tank thermal insulation
- spray system set into side for easy access to flask baskets
- integrated hand spray and spray protector
- all stainless steel construction
- visual and audible water-level indicator
- easy to clean without tools
- durable and energy-saving

2 Safety Symbols used in this Manual



Warning!

This is a warning of risk situations and dangers.

Failure to observe this warning could be life-threatening. These warnings has to be observed.



Information!

This symbol draws your attention to specific features that has to be observed.

3 Safety Guidelines

Configuring and operating this equipment requires ***precise knowledge*** and ***observance*** of the instructions in this ***user manual***. The equipment is designed only for its intended application.



WARNING:

Servicing and repairs should be carried out only by authorised specialists.
Disconnect the power plug before starting any maintenance work.



Make sure that the equipment is connected to the correct power source.



Attention! Danger of scalding!

Do not touch the hand shower at the nozzle or hose.

Do not point the spraying jet towards humans or animals.

The water, which comes out of the nozzle can reach temperatures of 95°C;
therefore only use the hand shower in the area of the screen cover or the hand basin.

4 Responsibility for Operation or Damage

The responsibility for operating the device lies exclusively with the owner or user if said device is incorrectly serviced, maintained or altered by persons not employed by an authorised dealer or if the device is used in a manner contrary to its specified purpose.

The unit has to be maintained and operated in accordance with this user manual.

Wassermann Dental-Maschinen GmbH is not responsible for damage arising from the **nonobservance** of these instructions.

Warranty and responsibility provisions contained in the sales and supply conditions of Wassermann Dental-Maschinen GmbH are not extended by these instructions.

5 Application

The Wapo-Ex 8 III is suitable for wax-scalding of up to 8 flask-halves as well as for scalding small parts on the low-lying perforated screen cover with the integrated hand shower.

Only use the device for this type of application.

6 Before Starting

6.1 Transport

Before transporting the unit, ensure that it has been unplugged from the power socket. Make sure that it is packed correctly in order to avoid accidental damage.



Be sure to check for any transport damage when unpacking the goods. Note down any damage if found.

6.2 Installation

Open the box, remove the packing materials, and carefully lift out the device and accessories. Check the included accessories.

The device has to stand horizontally on a steady and even surface



Install the device in a place where it will not block the working area and the functionality (take the dimensions into account).



Do not install the unit outdoors or in places without proper ventilation.

Before start-up, be sure the device reaches room temperature.



We strongly recommend to install the unit firmly to the discharge system, since it is equipped with an automatic water overflow discharge system, i.e. too much filled water runs out automatically, even when the discharge cock is closed. Furthermore complete emptying for cleaning purposes etc. is thus simplified.

6.3 Storage

The unit location is very important when it comes to workplace safety, even if it is only to be set up there for a short period. The room should be dry, well ventilated and vibration-free. An even temperature and wooden supports also help. The unit should not be stored or set up outdoors. If the unit is to be stored for an extended period, protect it from moisture and dust.


7 Installation/ Commissioning/ Menu Functions



Before commissioning, please install the following connections:


Firmly clamp the supplied drainage hose to the device and ensure there is a suitable drainage facility.

- Insert the plug into the socket – make sure that the voltage information between the device and the mains network is identical.
- The filter mat is already placed with the white side facing down.
- Fill the unit with water and **Wachs-Ex** up to the overflow limit (18.5 litres max.)
Dosage: approx. 40 ml / 18,5 litres water.
- Switch on the mains switch (green rocker switch).
- Symbol (A) must not appear.
- Also avoid an overflow!




Make sure you define the device's main parameters before switching it on. These settings are performed in the set-up menu (3) Familiarise yourself with the many options which the system offers before starting wax extraction. This applies in particular to functions (7) and (14)

7.1 The Set-Up Menu




7.2 The Set-Up Menu Page 1

Make the following settings in the set-up menu:




7.3 The Set-Up Menu Page 2





7.4 Description of Functions 5 - 14

5 Manual Wax Extraction Menu



The manual wax extraction process is initiated by touching the **wax extraction symbol**.


-  If the temperature falls below 63 °C the process is aborted and the heating symbol appears . The process must be manually started in order to continue working.

6 Automatic Wax Extraction Menu

Defines the parameters for the temperature and duration of the wax extraction process. By tapping on the numbers, the parameter selected is enlarged and the values can then be changed.



Touch the arrow to start the program. Interrupt or cancel the program using Start/Stop.



7 Defines the target temperature once the main switch has been activated.

Defines the target temperature once the main switch has been activated.

The pre-heat temperature also applies to waiting times and rest periods. Default: 65°C




If the programme temperature selected is higher, the temperature falls back to this value when the programme ends.



After defining the desired mode it may make sense **14** to adjust the pre-heat temperature of the programme temperature. This avoids long heating up or cooling down phases. The pre-heat temperature may be up to 8°C higher than the programme temperature for the programme to start. The control system will then adjust (reduce) the temperature to the programmed value.


8 Timer

This defines the weekday activation and deactivation times.



Programs running currently have a higher priority; they are not switched off until complete.

10 Display of the sleep function



Press the Sleep button to switch the system to sleep mode.

⑪ Definition of Calibration Method





DEF: Default (calibrated for mains water)


CAL Calibration (only necessary for extreme deviations in water quality)




This measure defines the conductance of the water!



⑫ Choose between °C Celsius and °F Fahrenheit




13 Pictogram for the Sequence in Time and Date Function



14 Heating ramp

Defines whether the time should start immediately or once the heating value has been reached.



(a) (Default) Unit heats up to the temperature programmed via 6 plus the tolerance to compensate for heat loss from the supply pump. After reaching this value, the result measured is validated by a control cycle. If the value is confirmed, the wax extraction process will start.

(b) The wax extraction process programmed begins immediately but not before reaching 65 °C Celsius. Please select a suitable preset temperature. 7

8 Operation

8.1 General Operating Instructions

All instructions for using the unit, whether in verbal or written form, are based on our own experience and experimentation and can only be regarded as guidelines.



The Wapo-Ex 8 III is protected against dry running in order to prevent the heating from burning out (flashing symbol). While the program is running, an alarm sounds when there is insufficient water. Take care that the water level is always above the white plastic ring (see in the housing), as the heating and the program will otherwise shut down automatically!

An additional safeguard protects the heating against overheating; it is automatically deactivated at above 97°C.



Wear heat-retardant gloves when working with the Wapo-Ex 8 III.




Unplug the device if it is not in use for a longer period.



Avoid overfilling and check the overflow system / drainage system. Strong overfilling can damage the electronics in the medium term. For damage caused by improper handling, no liability is assumed.

8.2 Operation

- Fill the unit with water and **Wachs-Ex** up to the overflow limit (18.5 litres max.)
Dosage: approx. 40 ml / 18,5 litres water.
- Also avoid an overflow!
- Activate main switch (green rocker switch).
- Check the water level – the low water alarm (A) must not flash.
- Place prepared cuvettes in the device clamp.
- Automatic wax extraction (temperature/time) is started by the (6) symbol. After the wax extraction process has ended, the temperature falls/rises to the preset temperature again. (7)
- For manual wax extraction, touch the (5) symbol. The process must be ended manually.



8.3 Hand Shower

In its ready state (> 65 °C) you can activate the hand shower at any time by moving the ball cock. The settings are as described in section 6.2. Always return the ball cock to its base position after use, i.e. in a horizontal position facing back. Please take note of the warning indications.



The Wapo-Ex 8 III is not ready to use until the temperature exceeds 65 °C, i.e. the pump will not start operating until the desired temperature is reached and also the count down of the timer do not start earlier. On addition there is a safety step: the heating stops automatically when the temperature of the heating is higher than 97°C.



Attention! Danger of Scalding!

Do not touch the hand shower at the nozzle or hose.

Do not point the spraying jet towards humans or animals.

The water, which comes out of the nozzle can reach temperatures of 95°C; therefore only use the hand shower in the area of the screen cover or the hand basin.



During scalding with the hand shower you must shift the ball cock (situated at the right side at the top of the housing) to the vertical position.

To relieve the hose and hand shower of pressure when not in use, the ball cock should always be returned to the horizontal position.

9 Troubleshooting

Fault	Cause	Solution
Wapo-Ex 8 III does not heat up	Dirty level sensor	Clean level sensor
	Water level too low	Add water
Water-level LED comes on and warning tone sounds. Fig. 1.	Dirty level sensor	Clean level sensor
	Water level too low	Top up water until control light goes off
	Re-calibrating the touch display	Turn the unit on and while the Wassermann logo appears tap the touch display twice. This will take you to the self-explanatory touch panel calibration process. Fig.: 2-3
Programme not starting although set temperature has been reached.	In mode (a) the unit intentionally overheats by 1°C and waits briefly to validate the value.	This is not a fault, please wait roughly 45 sec.



If the above recommendations do not solve the problem, contact your dental depot or our service department.




Fig. 1



Fig. 2




Fig. 3

10 Maintenance

10.1 Cleaning



Disconnect the power plug before starting any maintenance work.

The identification plate has always to be kept in easily legible condition and has not to be removed.



Remove external dirt from time to time with some form of cold cleaner.

Use only cold cleaners to avoid damaging the paintwork or removing the lettering.

The equipment should be cleaned at regular intervals to ensure trouble-free operation.

Once a week

- Change water and add **Wachs-Ex**.
- Dosage: approx. 40 ml / 18,5 litres water
- Clean the level tracer (white plastic ring), i.e. with a nail-brush.

Every 2 weeks

- - Replace the filter mat (white side facing down).

Every 4 weeks

- clean the unit with **Kalk-Ex**.
- Dosage: approx. 300ccm Kalk-Ex per full tank (18.5 litres water).
- If there is considerable limescale present, increase the amount of Kalk-Ex appropriately.
- The *initial supply* of Kalk-Ex comes with the machine.



Warning:

Leave Kalk-Ex to work for a maximum 12 hours.

Then drain the tank and rinse thoroughly.



Warning!

Please clean the water level sensor at least every 14 days with a brush.



10.2 Maintenance

The Wapo-Ex 8 III does not require any servicing. Just make sure that the device is kept clean.

10.3 Repairs



Servicing or repairs to the unit has only to be carried out by qualified technicians. Only original spare parts are to be used. Responsibility for the product is voided if unauthorised persons alter it or if inappropriate components are installed.

10.4 Spare Parts

If necessary please contact our service hotline phone.


10.5 Service Hotline 0049 (0)40 730 926 -20/ -24

10.6 Scope of Delivery/ Accessories

Included parts	Item no.:
Service set (4 filters, 1 l lime remover, 2 l wax remover)	170021
Discharge hose, PVC hose 19 x 4 (3/4"), 1.5 m	320030
Hose clip Ø 20 - 32 x 12.2	322017

Accessories:	Item no.:
Service set (4 filters, 1 l lime remover, 2 l wax remover)	170021
Mobile stand	170120
Discharge hose, PVC hose 19 x 4 (3/4"), per meter	320030

10.7 Warranty

 The warranty period for our equipment is 12 months. If faults occur within the warranty period, contact your dental depot or get in touch directly with our service department.

Your equipment should only be operated in perfect condition. If faults occur which could harm operators or third parties, the unit should not be used until it has been fixed.


This warranty does not cover damage caused by improper use, external mechanical causes, transport damage or interference with the unit by unauthorized persons.

11 Technical Data

Wapo-Ex 8 III	Item no.: 170981
Voltage*:	220–240 V / 50/60 Hz
Power consumption:	9 A
Output:	2000 W
Overall dimensions: Width x Height x Depth	610 x 540 x 420 mm
Weight:	42,0 kg
Noise emission:	≤ 70 dB (A)
Capacity:	18,5 l

*Other voltages on request.

The noise level of the unit amounts to ≤ 70 dB (A).

 Technical changes reserved.

12 Disposing of the Unit

The unit has to be disposed by an authorized recycling operation. The selected company has to be informed of all possibly health-hazardous residues in the unit.

12.1 Information on Disposal for Countries within the EU



To conserve and protect the environment, prevent environmental pollution and improve the recycling of raw materials, the European Commission adopted a directive that requires the manufacturer to accept the return of electrical and electronic units for proper disposal or recycling.

Within the European Union units with this symbol should not therefore be disposed of in unsorted domestic waste.

For more information regarding proper disposal please apply at your local authority.

13 EU Declaration of Conformity

in accordance with 2014/35/EU (Low Voltage Directive) and 2014/30/EU (Electromagnetic Compatibility Directive) and 2006/42/EC (Machinery Directive) and 2011/65/EU (RoHS Directive)

Manufacturer: **W A S S E R M A N N** Product description: Boil Out Unit for dental application
Dental-Maschinen GmbH
Rudorffweg 15-17
21031 Hamburg
Germany

Machine type: Wapo-Ex-8 III Item no. 170981

Applicable standards
DIN EN 61010-1
DIN EN 61000-6-3
DIN EN 61000-6-1
DIN 45635-1
DIN EN 60335-1

Hiermit wird bestätigt, dass die oben bezeichnete Maschine den genannten EU-Richtlinien entspricht. Diese Erklärung wird ungültig, falls die Maschine ohne unsere Zustimmung verändert wird.

This is to confirm that the above mentioned machine complies with the described EU guidelines. This declaration becomes invalid if the machine is modified without our approval.

Cette machine est conforme aux normes en vigueur d'Union européenne. Cet avis est nul et non avenant si cette machine est modifiée sans notre accord.

Esta máquina, anteriormente mencionada, cumple con los límites requeridos por el reglamento UE. Ahora bien, esta declaración quedará invalidada en caso de realizar modificaciones al aparato sin nuestra aprobación.

Hiermee wordt bevestigd dat bovengenoemde machine voldoet aan de voorgeschreven EU normen. Deze verklaring verliest geldigheid als er zonder onze uitdrukkelijke toestemming wijzigingen aan de machine worden aangebracht.

Place, date: Hamburg, 2017-07-06

Company stamp:



Signature: _____
Wilfried Wassermann
(Managing Director)



Wassermann Dental-Maschinen GmbH
Rudorffweg 15-17 · 21031 Hamburg · Germany
Phone: +49 (0)40 730 926-0 · Fax: +49 (0)40 730 37 24
info@wassermann-dental.com · www.wassermann.hamburg