WASSERMANN
Hamburg · Germany

Schleifscheibenwechsel

HSS-88

Diamantschleifscheibe

Replacement of the Grinding Disc

HSS-88

Diamond Grinding Disc

Changement de la meule

HSS-88

La meule diamantée

Hamburg · Germany
WASSERMANN
Since 1927

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GB USER MANUAL

FR MODE D´EMPLOI

DE Firmenanschrift:

Wassermann Dental-Maschinen GmbH
Rudorffweg 15-17
21031 Hamburg, Deutschland
Tel.: +49 (0)40/ 730 926 - 0 Fax.: +49 (0)40/ 730 37 24
sales@wassermann-dental.com
service@wassermann-dental.com
www.wassermann.hamburg
www.facebook.com/WassermannDentalmaschinen
Service-Hotline 040/ 730 92 6-50

GB Company address:

Wassermann Dental-Maschinen GmbH
Rudorffweg 15-17
21031 Hamburg, Germany
Phone : +49 (0)40 730 926 -20/ -24 Fax.: +49 (0)40 730 37 24
export@wassermann-dental.com
www.wassermann.hamburg
www.facebook.com/WassermannDentalmaschinen
Service Hotline 0049 (0)40 730 926 -20/ -24

FR Adresse :

Wassermann Dental-Maschinen GmbH
Rudorffweg 15-17
21031 Hambourg, Allemagne
Tél. : +49 (0)40 730 926 -20/ -24 Fax. : +49 (0)40 730 37 24
export@wassermann-dental.com
www.wassermann.hamburg
www.facebook.com/WassermannDentalmaschinen
Service d'aide : 0049 (0)40 730 92 6-24

.

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

GB Area of application

A) This guide describes how to replace an older diamond disc, a corundum disc, or a "Dia Quick" disc (Figure 2) with the current diamond disc (Figures 3 & 4).


B) This guide describes how to replace a current diamond disc with another current diamond disc (Figures 3 & 4).

Figures 2: older Disk

Figures 3: System 2026

Figures 4: System 2026

GB A) Replacement of the Diamond Grinding Disc HSS-88 (old > new)

Disconnect the power plug before starting any maintenance work.

1. Remove the front panel by unscrewing the four star-shaped knobs on the back. (If necessary, wrap a cloth around the star-shaped knob and use water pump pliers.)
2. Loosen the screw (4) using a size 13 ring spanner. To do this, wrap your right hand round the edge of the trimmer and press the cutting disc backwards with the palm of your hand.
3. If necessary, loosen the pressure pad (Z) by inserting a screwdriver between it and the grinding disc, and remove it together with the PVC disc (Y).
4. If necessary, clean the thread that is now visible using a wire brush.
5. If necessary, spray on some rust remover and leave it to work. tensioning element must be strong fitted).**PLEASE NOTE: Make sure you screw it in straight!**
6. Please carefully unscrew the puller supplied with the replacement grinding disc as far as possible. CAUTION: Make sure you screw it in straight!
7. Then, by turning the central screw, pull the grinding disc off the spindle. Do not use force; if necessary, spray on some more rust remover and leave it to take effect. To do this, tilt the trimmer backwards so that the grinding disc is facing upwards.
8. If necessary, polish the motor shaft to a shine using fine-grit sandpaper.
9. Pry the brass bushing (X) loose using two large screwdrivers. The brass bushing (X) and the pressure pad (Z) und the PVC disk (Y) are no longer required.
10. Remove any gypsum deposits from the bearing shell.
11. Lightly grease the bore of the new grinding wheel from the rear using the sealing grease supplied.
12. Push the new diamond disc onto the motor shaft as far as it will go.
13. Fit the clamping element without applying **grease.**
 (Note: Installation position! See figure 1 on page 4)
14. Screw the pressure bush (3) (lightly grease the thread) into the diamond disc and tighten firmly (secure the clamping element!). To do this, grip the edge of the trimmer with your left hand and press the grinding disc backwards.
15. Screw the bolt (4) and washer (5) into the end of the shaft in the same way and tighten firmly.
16. Check the grinding wheel for runout. To do this, look down at the grinding wheel from above and slowly turn it with your fingers. Normally, the grinding wheel runs smoothly! If the surface gives way at any point, use a sturdy tool to lever it back into position. Check the runout again until the wheel runs evenly.
17. Clean the sealing surfaces of the front panel and bearing shell and apply the new grease supplied.
18. Fit the front panel and tighten the lower star-shaped handles first, then the upper ones!
19. Remove any excess fat.

Please also use the delivered grease for greasing the windings!

GB B) Replacement of the Diamond Grinding Disc HSS-88

Disconnect the power plug before starting any maintenance work.

1. Remove the front panel by unscrewing the four star-shaped knobs on the back. (If necessary, wrap a cloth around the star-shaped knob and use water pump pliers.)
2. Loosen the screw using a size 13 ring spanner. To do this, wrap your right hand round the edge of the trimmer and press the cutting disc backwards with the palm of your hand.
3. Unscrew the pressure bushing.
4. Mount the puller for removing the old disc.
5. Grease the motor shaft (***not*** the first 25 mm from the front, because the grip-spring tensioning element (1) must be strong fitted).
6. Push the new diamond disc on the motor shaft up to the fixing.
7. Mount the grip-spring tensioning element (1) ***without*** grease. **(Attention: fitting position must like Figures 1, Site 4)**
8. Thread the pressure bushing (3) (grease the winding a little bit) into the diamond disc and fix (**clamp the grip-spring tensioning element! (1)**).
9. Thread the screw (4) (grease the winding) with disc (5) into the and of the shaft and screw up.
10. Screw the pressure bush (3) (lightly grease the thread) into the diamond disc and tighten firmly (secure the clamping element!). To do this, grip the edge of the trimmer with your left hand and press the grinding disc backwards.
11. Check the grinding wheel for runout. To do this, look down at the grinding wheel from above and slowly turn it with your fingers. Normally, the grinding wheel runs smoothly! If the surface gives way at any point, use a sturdy tool to lever it back into position. Check the runout again until the wheel runs evenly.
12. Clean the surfaces of the front-plate and bearing box and grease with the delivered grease. Do not spread it.
13. Put back the front plate and fasten the star handles tightly. First the bottom star handles, then the top star handles. Now spread the grease and remove any grease.
14. Remove any excess fat.

Please also use the delivered grease for greasing the windings!

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