

**Fine Blasting Units Wassermann Dental-Maschinen GmbH**  
**Cemat-4 II, -3 II, -2 II/ Gobi Comfort-4, -3, -2/ Gobi Classic/ Puna II**

**Nozzles and Blasting Material**

Colour code* (tank/ handpiece):	yellow	blue	green	red
Fine blasting nozzle (handpiece)	0,6	0,8	0,8	1,2
Nozzle (tank) 1,4 ( < 110 µm)	X	X		
Nozzle (tank) 2,0 ( > 110 µm)			X	X
Aluminium Oxyde 50 µm	+++	++	+	-
Aluminium Oxyde 110 µm	+	+++	+++	+
Aluminium Oxyde 250 µm	-	-	+	+++
Glass Beads 50 µm	++	+	-	-
Glass Beads 110 µm	+	++	+++	+
Glass Beads 125 µm	-	+	++	+++

**Annotation:**

- X = Combination fine blasting nozzle (handpiece)/ nozzle (tank)
- = does not fit
- + = fits
- ++ = fits well
- +++ = optimal

The ratings in the table are standard values that may vary depending on the nature of the blasting material.

**\*This colour code with the above mentioned combination of the nozzles is delivered as standard:**

Item no. 190984 Cemat-4 II	Item no. 190884 Gobi Comfort-4	Item no. 190881 Gobi Classic	Item no. 190885 Puna II
Item no. 190983 Cemat-3 II	Item no. 190883 Gobi Comfort-3		
Item no. 190982 Cemat-2 II	Item no. 190882 Gobi Comfort-2		
Cemat-4 II incl. 4 fine blasting nozzles: 0,6 (yellow), 0,8 (blue), 1,2 (red), 0,8 (green)	Gobi Comfort-4 incl. 4 fine blasting nozzles: 0,6 (yellow), 0,8 (blue), 1,2 (red), 0,8 (green)		
Cemat-3 II incl. 3 fine blasting nozzles: 0,8 (blue), 1,2 (red), 0,8 (green)	Gobi Comfort-3 incl. 3 fine blasting nozzles: 0,8 (blue), 1,2 (red), 0,8 (green)		
Cemat-2 II incl. 2 fine blasting nozzles: 0,8 (blue), 1,2 (red)	Gobi Comfort-2 incl. 2 Fine blasting nozzles: 0,8 (blue), 1,2 (red)	Gobi Classic incl. 2 fine blasting nozzles: 0,8 (blue), 1,2 (red)	Puna II incl. 2 fine blasting nozzles: 0,8 (blue), 1,2 (red)