

HSS-88 incl. grinding discs	Item no.
HSS-88 incl. diamond grinding disc	168998
HSS-88 incl. Dia-Quick grinding disc	168996
HSS-88 incl. corundum grinding disc	168999
HSS-88 incl. foil grinding disc	168991



HSS-88 with inclined support, removable



Extendable by orthodontic model table and/or laser light guide for 3 D trimming, ideal for orthodontics (the orthodontic model table can be delivered also in Swiss and EOS version)

## MADE IN GERMANY

Wassermann is a medium-sized, owner-operated company with head office in Hamburg. We design and produce high-quality technical devices in Germany. Our main users are dental laboratories, dental clinics and various industrial areas worldwide.

### Why a Wassermann product?

Wassermann products are manufactured from high-quality materials and the proportion of plastics used is intentionally kept as low as possible. Our devices are known for their suitability for daily use, functionality and longevity. As a result, the users place high levels of trust in the Wassermann brand. We feel it is our duty to support our customers on a daily basis when it comes to safety, service and productivity. Technological perfection, optimum quality and especially sustainability are the benchmarks we set. Thanks to the combination of these values, the family-owned company, since its establishment in 1927, has gained worldwide recognition as a manufacturer of high-quality products.



Further information are available on our website at [www.wassermann.hamburg](http://www.wassermann.hamburg)



/WassermannDentalmaschinen

Your Dental Dealer

## HSS-88



## Trimmer



Wassermann Dental-Maschinen GmbH  
Rudorffweg 15-17 · 21031 Hamburg · Germany  
Phone: +49 (0)40 730 926-20/ -24 · Fax: +49 (0)40 730 37 24  
export@wassermann-dental.com · www.wassermann.hamburg  
www.facebook.com/WassermannDentalmaschinen

# Trimmer HSS-88

The classic among wet trimmers - powerful, flexible and durable: The multi-faceted HSS-88 by Wassermann is the most sold trimmer for good reason. Its high torque and high speed, but also energy efficient motor guarantees the highest grinding speed. This maintenance-free motor does not require fans or cooling fins, since the motor is cooled through the outer casing. And of course it is needless to say that the trimmer operates with very little noise and no vibration.



## Major Features

- Very powerful, maintenance-free motor, max. 3000 min<sup>-1</sup>, splash protected according to IP 54
- Sturdy, oxidation-free bench with longitudinal slots and side openings – avoids water spills
- Supplied with solenoid valve for automatic water supply
- Multiple expansion options
- **Included in delivery:**  
Inclined support  
Supply hose, 1 m  
Discharge hose, 1 m

Accessories	Item no.
Spare diamond grinding disc, incl. puller and grease	168063
Spare Dia-Quick grinding disc (partially diamond coated), incl. grease	162063
Spare corundum grinding disc, incl. puller and grease	168024
Spare foil grinding disc, 1 set = 4 pieces, incl. grease	162013
Spare foil grinding disc, fine, 1 set = 4 pieces, incl. grease	162023
Conversion kit, Dia-Quick and foil grinding disc	168071
Pressing tool to roughen corundum grinding disc	710001
Sealing grease for trimmer, 25 ml	629002
Puller for grinding disc change	164026
<b>Orthodontic accessories for producing 3D orthodontic models</b>	
Orthodontic model table KMS according to Begg	168140
Laser light guide for model table, cpl.	168120
Orthodontic template for the upper jaw (to grind the lower jaw/tuber level)	164034
<b>HSS-88 without grinding disc</b>	
	<b>Item no. 168997</b>
Voltage	220–240 V / 50/60 Hz
Power consumption	4.6 A
Output	1000 W
W x H x D	330 x 335 x 410 mm
Weight	19.8 kg
Sound level	≤ 70 dB (A)
Speed	3000 min <sup>-1</sup>
Duty type	S6-40 %
Table dimensions projecting	120 mm