



MADE IN GERMANY

Wassermann is a medium-sized, owner-operated company with head office in Hamburg. We design and produce high-quality technical devices in Germany. Our main users are dental laboratories, dental clinics and various industrial areas worldwide.

Why a Wassermann product?

Wassermann products are manufactured from high-quality materials and the proportion of plastics used is intentionally kept as low as possible. Our devices are known for their suitability for daily use, functionality and longevity. As a result, the users place high levels of trust in the Wassermann brand. We feel it is our duty to support our customers on a daily basis when it comes to safety, service and productivity. Technological perfection, optimum quality and especially sustainability are the benchmarks we set. Thanks to the combination of these values, the family-owned company, since its establishment in 1927, has gained worldwide recognition as a manufacturer of high-quality products.



Further information are available on our website at www.wassermann.hamburg



/WassermannDentalmaschinen

Your Dental Dealer _____

Wapo-Mat III



Polymerisation Unit



Wassermann Dental-Maschinen GmbH
Rudorffweg 15-17 · 21031 Hamburg · Germany
Phone: +49 (0)40 730 926-20/ -24 · Fax: +49 (0)40 730 37 24
export@wassermann-dental.com · www.wassermann.hamburg
www.facebook.com/WassermannDentalmaschinen

Polymerisation Unit

Wapo-Mat III

The water bath of the **Wapo-Mat III** can be used for long-time and boiling polymerisation of all common plastics. The user set the relevant parameters quickly and easily using self-explanatory symbols. The option of 3 temperatures and 3 times is possible, with remaining time reading. The factory setting is that the main polymerization process will start automatically at 95 °C ninety minutes before the program ends.

Process up to
9 flasks in a flask holder
at one time

Energy saving
thermal insulation



Touch screen for
fast and easy operation

Major Features

- Compact benchtop automat with fully automatic program sequence
- Powerful, electronically controlled heater for optimal water temperature and long service life of the electrical components
- Low water level indicated visually and audibly
- Low maintenance, easy to clean
- **Included in delivery:**
 - 1 L lime remover
 - Discharge hose, 1.5 m
 - Hose clip Ø 20-32 x 12.2

Accessories	Item no.
1 L lime remover	655002
Discharge hose per meter	320030

Wapo-Mat III	Item no. 170980
Voltage	220-240 V / 50/60 Hz
Power consumption	8 A
Output	1800 W
W x H x D	330 x 460 x 530 mm
Weight	18.5 kg
Sound level	≤ 70 dB (A)
Water capacity	18.5 l