



Steam gun



Splash protection made of stainless steel with plastic basement, available as accessory

MADE IN GERMANY

Wassermann is a medium-sized, owner-operated company with head office in Hamburg. We design and produce high-quality technical devices in Germany. Our main users are dental laboratories, dental clinics and various industrial areas worldwide.

Why a Wassermann product?

Wassermann products are manufactured from high-quality materials and the proportion of plastics used is intentionally kept as low as possible. Our devices are known for their suitability for daily use, functionality and longevity. As a result, the users place high levels of trust in the Wassermann brand. We feel it is our duty to support our customers on a daily basis when it comes to safety, service and productivity. Technological perfection, optimum quality and especially sustainability are the benchmarks we set. Thanks to the combination of these values, the family-owned company, since its establishment in 1927, has gained worldwide recognition as a manufacturer of high-quality products.



Further information are available on our website at www.wassermann.hamburg

Your Dental Dealer

Wasi-Steam Classic



Steam Cleaner



Wassermann Dental-Maschinen GmbH
Rudorffweg 15-17 · 21031 Hamburg · Germany
Phone: +49 (0)40 730 926-0 · Fax: +49 (0)40 730 37 24
info@wassermann-dental.com · www.wassermann.hamburg

Steam Cleaner

Wasi-Steam Classic

Compact format, strong performance: The **Wasi-Steam classic** is the perfect solution when there is little space and big demand. With its variable work pressure of up to 4.2 bar it offers in particular the multi-station user a gentle but quick and thorough cleaning of all small parts.



Manual filling,
flexible choice of location

Comfortable usage because of
single-handed operation
of the steam gun

Variable working pressure

Thermo-safety switch
with manual reset

Major Features

- Perfect for use with limited space
- Adjustable up to 4.2 bar
- Short heat-up phase because of 1,200 Watt-special steel heating
- Gentle and quick to use, thorough results
- High safety because of protection against dry running
- Wall mounting possible
- Tank and diverse housing parts made of stainless steel
- **Included in delivery:**
Fastening set for wall mounting

Accessories	Item no.
Splash protection (W 235 x H 385 x D 294 mm)	174360
1 L lime remover	655002

Wasi-Steam Classic	Item no. 174993
Voltage	220–240 V / 50 / 60 Hz
Power consumption	5.3 A (6.3 A medium-blowing-fuse)
Output	1220 W
W x H x D	235 x 340 x 340 mm
Weight	11.0 kg
Sound level	≤ 70 dB (A)
Working pressure	max. 4.2 bar
Steam temperature	approx. 151°C
Boiler capacity	3.3 l
Safety class	IP 22