

Space-saving & powerful –
the perfect combination.

easy · powerful · versatile

Wasi-Steam 3 delivers reliable steam power with consistently high performance. Ideal for cleaning, scalding, and degreasing in the lab.

Benefits: Fits the existing wall bracket from Wasi-Steam 2 – no new drilling needed. Maintains stable pressure at 6 bar and 164°C. Durable stainless steel housing designed for continuous operation. **New:** LED indicators for all operating states. Enhanced safety with automatic monitoring of water supply, leakage, and low water levels.

OUR PROVEN SUCCESS FORMULA:

We combine innovation, expertise, and teamwork.



Always here for you

Need a quote, want to place an order, or have a question about delivery, service, or warranty?

We're happy to take your call and help you personally.



www.wassermann.hamburg



fb.com/WassermannDentalmaschinen

Steam Cleaner Wasi-Steam 3



Fig.: Wasi-Steam 3
+ Splash Guard Unit

Info-Hotline

+49 (0)40 730 926-0



SCAN ME 



Fig.: Wasi-Steam 3
+ Splash Guard Unit

WASSERMANN
Hamburg · Germany

Wasi-Steam 3 Steam Cleaner



QUALITY MADE IN GERMANY SINCE 1927

The professional steam cleaner:

High-performance steam unit with fixed water connection,
ideal for all lab cleaning tasks. Robust stainless steel
construction – durable, sustainable, and hygienic.
Easy to use with a high level of safety.
Flexible use as a benchtop or wall-mounted unit.



Key features at a glance:

- Selectable dry or wet steam
- LED indicators for all operating states
- Compact design for benchtop or wall mounting
- Steam pressure up to 6 bar at 164°C
- Automatic water refill via integrated pump
- Standby mode reduces power use and protects components
- Low maintenance and long service life

Optional: steam quantity control and manual refill system

Included in delivery:

Pre-installed supply and drain hoses,
wall-mounting kit, steam gun holder

Technical data:

Power: 1,270 W
Current: 5.5 A
Tank capacity: 3.3 l

Dimensions (W × H × D max.): 505 × 400 × 230 mm
(+ Steam gun hose)
Housing dimensions (W × H × D): 400 × 400 × 230 mm

Steam pressure: 6 bar
Weight: 23 kg
Sound level: ≤ 70 dB(A)

Wasi-Steam 3
Item no.: 174989

Wasi-Clean Water Softening Unit
Item no.: 174995

Splash Guard Unit, stainless steel
Item no.: 174230

Wall-Mount Holder for Splash Guard
Item no.: 174625

Special Nozzle D 2.5 mm
Item no.: 174352

Refill Unit (5 liters)
Item no.: 174430

Steam Quantity Regulator Kit
Item no.: 174633

Factory-installed
option available

Info & Support:

Wassermann Dental-Maschinen GmbH
Rudorffweg 15 - 17
21031 Hamburg, Germany

Tel: +49 (0)40 730 926-0
www.wassermann-dental.com
sales@wassermann-dental.com