

# The new way to capture grinding dust compact, mobile, powerful

## Airbox One

combines a dust box and extraction unit in one – for clean, safe working directly in the treatment room and the lab.

**Advantages:** An end to airborne grinding dust in the dental practice. Twin-fan configuration for maximum suction, plus powerful LED lighting for a clear view of your work.

**New:** Mobile extraction with dust class M filtration and particle removal efficiency of up to 99.99%.

### OUR PROVEN FORMULA FOR SUCCESS:

We combine technology, experience, expertise and team spirit.



### Always here for you

Would you like to request a quote, learn more about our products, or do you have questions about your order, delivery or warranty?

We look forward to your call and will be happy to help.



[www.wassermann.hamburg](http://www.wassermann.hamburg)



[fb.com/WassermannDentalmaschinen](https://www.facebook.com/WassermannDentalmaschinen)

## Airbox One Dust Extraction Box



Fig.: Airbox One with protective covers

### Info Hotline

+49 (0)40 730 926-0



SCAN ME



Fig.: Airbox One with protective covers

# Airbox One

## Dust Extraction Box



**QUALITY MADE IN GERMANY SINCE 1927**

### The compact professional extraction system:

Dust box and extraction unit **combined in one device.**

Flexible portable use from the treatment room to the dental laboratory.

Clean and safe thanks to dust class M filtration.

Perfect for short grinding jobs and and portable use.

### The key features at a glance:

- Dust box and extraction unit combined in one unit
- Powerful and quiet
- Practical carrying handle
- Class M filter with a particle removal efficiency of up to 99.99%
- Simple LED lighting for an optimal view of your work
- Full mobility for use in any room
- Ergonomic shape for optimal posture
- Tilt and slip resistant, durable and easy to clean
- High-quality finish, made of stainless steel

**Optional:** Protective shields with pass-through sleeves and separate openings for handpieces, or arm supports

### Technical data

Power:	65 W	Voltage:	100-240 V / 50/60 Hz	Noise level:	60 dB
Current:	1 A	Dimensions W x H x D:	370 x 320 x 370 mm	Max. airflow:	1200 l/min
		Weight:	6.8 kg	Max. pressure:	700 Pa

#### Airbox One

incl. daylight LED lighting  
Item no.: 143372

#### Set of protective covers

Item no.: 143305

#### Magnifier

Item no.: 143251

#### Set of arm supports

Item no.: 143303

### Info & Advice:

Wassermann Dental-Maschinen GmbH  
Rudorffweg 15 - 17  
21031 Hamburg, Germany

Tel: +49 (0)40 730 926-0  
www.wassermann-dental.com  
export@wassermann-dental.com